

Sanitize Your Computer Mouse

Supplies Needed to Clean and Sanitize Your Computer Mouse

- Isopropyl alcohol, 70%
- Q-Tips
- Toothpicks
- Microfiber Cloth



1. Unplug mouse from computer.
2. Wipe mouse down with dry cloth.



3. Run a toothpick around the crevices in the mouse's body and around the bottom of the mouse feet and sensor.



4. Dip your Q-tip in isopropyl alcohol and wipe down the dirty parts of your mouse including the crevices and sensors.



5. Allow mouse to dry completely before plugging back in.